

# Bias Free Linear Vector for Cloning Recalcitrant DNA & Accelerating Sequence Finishing

Ronald Godiska<sup>1</sup>, Rebecca Hochstein<sup>1</sup>, Sarah Vande Zande<sup>1</sup>, Nikolai Ravin<sup>2</sup>, Attila Karis<sup>3</sup>, David A. Mead<sup>1</sup>

<sup>1</sup> Lucigen Corporation, Middleton, WI 53562

<sup>2</sup> Centre Bioengineering, Russian Academy of Science, Moscow, Russia

<sup>3</sup> Mississippi State University


## Abstract

We have developed a novel linear vector for unbiased cloning of 0-30 kb inserts in *E. coli*. This vector, termed "pJAZZ", shows unprecedented ability to maintain large inserts from very AT-rich genomes. The otherwise difficult-to-clone genome from *Flavobacterium columnare* (70% AT, 3.2 Mb) was sequenced to seven fold coverage using the pJAZZ vector, with only 10 sequencing gaps. The linear vector was able to maintain 20-30 kb fragments from *Lactobacillus helveticus* (65% AT) and 2-4 kb inserts from *Piromyces* (up to 96% AT), which were unclonable in conventional plasmids. Unlike fosmid cloning, the construction of large-insert libraries (10-20 kb) in pJAZZ is simple and robust, using standard methods of transformation and plasmid purification. We are evaluating the use of a single pJAZZ shotgun library to eliminate the need for multiple libraries, making finishing easier and more cost effective. Enhanced stability of inserts in the pJAZZ vector is attributed to both the lack of supercoiling and the lack of transcriptional interference. Torsional strain inherent to supercoiled plasmids can induce localized melting and generate secondary structures, which are substrates for deletion or rearrangement by resolvases and replication enzymes. For example, the instability of tandem repeats and palindromic sequences is presumably due to cleavage of hairpin structures or to replication slippage across the secondary structures. Most conventional plasmid vectors also induce strong transcription and translation of inserted fragments, and they allow transcription from cloned promoters to interfere with plasmid stability. As a result, certain DNA sequences are deleterious or highly unstable, leading to sequence "stacking", clone gaps, or a complete inability to construct libraries, especially from AT-rich genomes or toxic cDNAs. The transcription-free, linear pJAZZ vector also minimizes "sequence gaps" caused by secondary structures, as shown by its stable cloning of inverted repeats and di- and tri-nucleotide repeats.

## Replication of the pJAZZ Linear Vector


The pJAZZ vector is derived from the linear phage N15 of *E. coli*. The *repA* and *telN* genes of N15 encode replicase and telomerase. Bi-directional replication is initiated by RepA and is carried out by the host DNA polymerase. TelN cleaves the replicated telomeres and covalently joins the 5' and 3' strands of each free end, re-creating the terminal hairpin loops (*telL* and *telR*).

Linear plasmid DNA is electroporated into BigEasy™ TSA *E. coli* cells, which contain N15 genes needed for partitioning and copy number regulation. Linear plasmid DNA is efficiently isolated using standard methods, such as alkaline lysis and binding to a silica matrix.




## Minimal Size Bias

Sheared *E. coli* genomic DNA was cloned into the pJAZZ vector or pUC19 without size selection. As expected, the size distribution of the pUC19 inserts was clearly skewed toward smaller inserts. In contrast, the distribution of the pJAZZ inserts closely matched the input DNA.



Size distribution of pJAZZ vs pUC19 clones. The size distribution of 85-100 inserts in each vector is graphed. The insert DNA and 1 kb ladder are shown below the chart.

## pJAZZ™ Linear Vectors for High Stability Cloning



The linear pJAZZ vectors lack the structural and lysis genes of N15. All versions of the pJAZZ vectors are also transcription-free, and they appear to behave similarly. Dual selection ensures that recombinant clones contain both arms of the vector. The origin of replication is used to select for the left arm in the pJAZZ-OC vectors.


The ends of the vectors are free to rotate during replication, so cloned inserts are not subject to torsional stress caused by supercoiling. Transcriptional terminators at the cloning site minimize

transcriptional interference between the insert and the vector, increasing insert stability. The pJAZZ vectors are low-copy (~5/cell) to further promote stable propagation of inserts, and their copy number can be induced 5-20X for DNA preparation.

The *E. coli* host strain for the pJAZZ vectors contains several required genes of N15: *telN* encodes telomerase for replication; *sopAB* is needed for stable maintenance, and *antA* regulates copy number.

## Cloning Large PCR products


PCR products of 15 - 20 kb were amplified from phage lambda DNA and cloned into the pJAZZ vector. Plasmid DNA from transformants was cut with *NotI* to excise the insert, then analyzed by gel electrophoresis. The left vector arm migrates at 10 kb; the right arm (2 kb) was run off the gel.







## Cloning "Unclonable" DNA

### AT-rich genomic DNA

Using conventional circular vectors, AT-rich DNA is often difficult to clone in *E. coli*, producing very few stable, intact clones (Figure A).




In contrast, large clones of extremely AT-rich DNA are stable in the linear, transcription-free pJAZZ vectors (Figures B, C, D, E)



DNA were sheared to 2-6 kb (Figures A & C) or 6-20 kb (Figures B, D, E), end-repaired, size-selected, and cloned into the vector. Plasmid DNA from pJAZZ transformants was cut with *NotI* to excise the insert and analyzed by gel electrophoresis. Vector bands are 12 kb and 2 kb. AT content of inserts was as high as 96% (Figure D).


### GC-rich genomic DNA

The GC-rich genome of *B. animalis* (67% GC) was sheared to 6-20 kb, end-repaired, and cloned into the pJAZZ vector. Clones from this library were similar in size and number to those from the AT-rich *L. helveticus* library (see below). Uncut DNA from transformants was 24-40 kb, corresponding to inserts of 10-26 kb.




## Repetitive DNA

Repetitive DNA sequences can be extremely difficult to clone in circular plasmids. For example, mollusk cDNA of 0.3-2 kb was unclonable in all circular plasmids tested. However, the pJAZZ linear vector produced a library of full length clones. Repetitive DNA from two of the clones is shown (Figures A & B).



In *E. coli*, clones of 50-100 CGG repeats from the Fragile X locus are highly unstable in circular vectors, and the frequency of deletion is increased by transcription and supercoiling. Corroborating these results, the pJAZZ vector was able to maintain fragments containing 220 copies of the CGG repeat, which has not been achieved with circular vectors.




## Efficient Library Construction and Assembly

The pJAZZ vector was used to construct genomic libraries of *Flavobacterium columnare* (70% AT, 3.2 Mb). Clones of 2-10 kb were readily obtained, and gave sequence reads of ~800 bp.

Automated assembly closely followed the theoretical Lander Waterman curve, yielding 21 major contigs. Further analysis revealed just 10 sequence gaps. Therefore, this AT-rich genome was completely cloned and nearly completely sequenced without the use of fosmid or BAC libraries.

-Attila Karsi, Mississippi State University



## Summary

- Linear, transcription-free pJAZZ vectors allow cloning of repetitive DNA and large AT-rich or GC-rich DNAs. Libraries of 10-20 kb clones of AT-rich genomic DNA were routinely created with the linear transcription-free vector. Likewise, repetitive DNAs were maintained without rearrangement.
- Efficient Genomic Cloning and Sequencing. Powerful, unbiased cloning greatly reduced the need for manual finishing of genomic libraries. A 3-Mb microbial genome was assembled without the use of fosmid clones.
- Rapid and simple protocol. No vector preparation or special techniques are needed to generate high quality, large-insert libraries of otherwise "unclonable" DNAs.

**Lucigen®** Advanced Products for Molecular Biology

Lucigen Corporation 2120 W. Greenview Dr., Ste 9 Middleton, WI 53562 www.lucigen.com