



# **Beyond Basic Target Enrichment:** **New Tools to Fuel Your NGS** **Research**

*Sequencing, Finishing and  
Analysis in the Future 2012*




**Agilent Technologies**


Not approved for use in diagnostic  
procedures

# Agenda


- **1** New tool for DNA methylation
- 2** Deep RNA sequencing
- 3** Fast strategy with benefits



# Differentially Methylated Regions (DMR)



- CpG islands
- CpG island **shores** (~2kb away from islands)
- CpG island **shelves** (~4kb away from islands)
- CpG shore methylation is strongly related to gene expression




HS3ST4 :heparan sulfate D-glucosaminyl 3-O- sulfotransferase 4

# Whole Genome Bisulfite Sequencing (WGBS)

- Whole Genome Coverage
- Costly and time consuming
- Requires extensive bioinformatics
- Limited scalability per run

## Target Enrichment:

- reduce costs/sample
- maintain similar whole-genome CpG coverage
- Have increased scalability





# SureSelect<sup>XT</sup> Human Methyl-Seq

## Discovery Tool

- Not methylation-state dependent
- No prior knowledge needed

## Comprehensive design

- Not limited to CpG Islands.
- Comprehensive Content
- CpG Islands, Promoters and DMRs


## DESIGN CONTENT - 84 Mb

### Design, 3.7M CpGs

- CpG islands
- Cancer, Tissue-specific DMRs
- GENCODE promoters
- DMRs or regulatory features in:
  - ✓ CpG Islands, shores and shelves  
±4kb
  - ✓ DNaseI hypersensitive sites
  - ✓ Refseq Genes
  - ✓ Ensembl Regulatory Features



# SureSelect Methyl-Seq Protocol



No PCR before  
Bisulfite treatment  
to preserve the  
Methylation state

Bisulfite treatment is  
performed after  
hybridization to  
maximize sample  
complexity

# Methyl-Seq Comparison with published WGBS data

Vol 462 | 19 November 2009 | doi:10.1038/nature08514

nature


## Human DNA methylomes at base resolution show widespread epigenomic differences

Ryan Lister<sup>1\*</sup>, Mattia Pelizzola<sup>1\*</sup>, Robert H. Downen<sup>1</sup>, R. David Hawkins<sup>2</sup>, Gary Hon<sup>2</sup>, Julian Tonti-Filippini<sup>4</sup>, Joseph R. Nery<sup>1</sup>, Leonard Lee<sup>2</sup>, Zhen Ye<sup>2</sup>, Que-Minh Ngo<sup>2</sup>, Lee Edsall<sup>2</sup>, Jessica Antosiewicz-Bourget<sup>5,6</sup>, Ron Stewart<sup>5,6</sup>, Victor Ruotti<sup>5,6</sup>, A. Harvey Millar<sup>4</sup>, James A. Thomson<sup>5,6,7,8</sup>, Bing Ren<sup>2,3</sup> & Joseph R. Ecker<sup>1</sup>

Cell line: IMR90 (female lung fibroblast)





# SureSelect vs. Whole-genome bisulfite sequencing



# Concordance with WGBS data

## Reproducibility






# Agenda

1 New tool for DNA methylation

2 Deep RNA sequencing

3 Fast strategy with benefits



# RNA-Seq Transcriptome Target Enrichment in Cancer

## SureSelect RNA



### Targeted next-generation sequencing of a cancer transcriptome enhances detection of sequence variants and novel fusion transcripts

Joshua Z Levin\*, Michael F Berger<sup>†</sup>, Xian Adiconis\*, Peter Rogov\*, Alexandre Melnikov\*, Timothy Fennell<sup>‡</sup>, Chad Nusbaum\*, Levi A Garraway<sup>†§</sup> and Andreas Gnirke\*


*Genome Biology* 2009, **10**:R115

□ **Aim:** RNA Target Enrichment of 467 Cancer Genes (~all Tyr Kinases + Genes from Cancer Gene Census); Overall >800 target transcripts

□ **Method:** Enrich cDNA from K-562 CML cell line cDNA libraries and compare results from before and after target enrichment



Agilent Technologies



Not approved for use in diagnostic procedures

# Benefits of RNA-Seq Target Enrichment in Cancer

## SureSelect RNA

Feature	Without RNA Target Enrichment	With RNA Target Enrichment	
Avg. Coverage	14.4x	606x	42X greater coverage
Fold-Enrichment	N/A	42x	
Regions with 20x coverage	13%	63%	
Common SNPs	76	257	More sensitive variant detection
Novel Variants	4	16	
Splice-Junctions	2958	4720	
Alternatively spliced-genes	52	177	
Gene Fusions	2	6	

- Increased coverage improves ability to find rare transcripts and novel fusions
- Whole transcriptome analysis requires >40x more sequencing to achieve same depth for targeted regions



Agilent Technologies




Not approved for use in diagnostic procedures

# Agenda

1 New tool for DNA methylation

2 Deep RNA sequencing

● 3 Fast strategy with benefits



Agilent Technologies

Not approved for use in diagnostic  
procedures

# HaloPlex – Simplicity of Next Generation PCR

HaloPlex

## SIMPLE WORKFLOW PREMIUM PERFORMANCE

- ✓ **Low input (200ng) DNA**
- ✓ **Library-prep free**
- ✓ **Less than 6 hours**
- ✓ **96 indexed samples for throughput**
- ✓ Simple Design Wizard
- ✓ Single tube prep
- ✓ Any Desktop Sequencer
- ✓ **Any Targets up to 5Mb**



### **Compared to SS:**

- **Faster to results due to no library prep and short**
- **Small panels**
- **Interrogates repeats**
- **Human only**

New!

New!



Agilent Technologies

Not approved for use in diagnostic procedures



# HaloPlex: 4-step protocol


HaloPlex

1. Digest DNA

2. Hybridize  
Probes



3. Purify / Ligate

4. Amplify



# HaloPlex with FFPE Samples and WGA Samples

HaloPlex



Cancer Panel developed by OICR  
19 Genes, ~60Kb size

Cancer Panel developed at Karolinska  
Institutet 78 Genes, ~197 Kb size



Agilent Technologies

Not approved for use in diagnostic  
procedures

# Complete NGS Workflow

Library Prep

Target  
Enrichment

Sequencing

Data Analysis



**SureSelect Kits**



**SureSelect Kits HaloPlex Kits**



**QC: Bioanalyzer/  
TapeStation**



**QC: Bioanalyzer/  
TapeStation**



**PCR: Surecycler 8800  
qPCR: Mx3005P**



**PCR: Surecycler 8800  
qPCR: Mx3005P**



**Agilent NGS Automation**



**Agilent NGS Automation**

**Compatible with:**



**Illumina GALLx/HiSeq**



**SOLiD/5500**




**Roche 454/GS Junior**



**Ion Torrent PGM/MiSeq**



**Pacific Biosciences**




**Agilent Technologies**

Not approved for use in diagnostic  
procedures

# SureSelect: Complete “Omics” Solution

- ✓ **DNA:** Genetic variation
- ✓ **RNA:** Gene Expression
- ✓ **Methylation:** Effects on Gene Expression





Thank you.



Agilent Technologies