



RapTOR: Automated sequencing library preparation and suppression for rapid pathogen characterization

Todd W. Lane, Principal Investigator

SFAF Meeting, June 7, 2012

Santa Fe, NM




Sandia is a multiprogram laboratory operated by Sandia Corporation, a Lockheed Martin Company, for the United States Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.





RapTOR system concept






RapTOR is designed to address large outbreak scenarios




Outbreak dynamics are measured in **days to weeks**.

Culture-based identification of agent can take **months to years**.




Probe-based methods are fast but often confounded.


- Don't recognize deeply diverged targets.
- Blind to unanticipated features (eg, IL-4).
- Unusual profiles are difficult to interpret.
 - Mix-and-match features? Artifact?




LLMDA



GreeneChip



ViroChip



Second Generation Sequencing methods are transformational for pathogen characterization



Brute-force SGS of clinical samples has enabled the discovery of novel pathogens




http://www.illumina.com/systems/hiseq_2000.ilmn

| Disease | Sample | Novel Agent Detected | Total Reads | Hits on Agent | Reference | Suppression for 1X coverage |
|-----------------------------------|----------------------|------------------------------|-------------|------------------------|--------------------------|-----------------------------|
| Merkel cell carcinoma | tumors | "Merkel cell polyomavirus" | 395,734 | 2 (0.00005%) | Science 319:1096 '08 | 14-18X |
| organ transplant related fatality | serum & organs | "Dandenong" arenavirus | 103,632 | 14 (0.014%) | N Engl J Med 358:991 '08 | 5X |
| pediatric gastroenteritis | feces | "human klassevirus " | 937,935 | 849 (0.09%) | Virol J 6:82 '09 | 0 (75%) |
| pediatric influenza-like illness | nasopharyngeal swabs | "human enterovirus type 109" | 20,825,810 | 119 (0.0006%) | J Virol 84:9047 '10 | 0 (40%) |


Deplete non-informative NA to improve efficiency of SGS analysis

Increasing pathogen signal to noise can be accomplished through suppression





The RapTOR approach was designed to address future challenges of SGS-based sample analysis



Sboner et al. *Genome Biology* 2011, 12:125

SFAF June 7, 2012

10 -100 X suppression would have significant effects on real samples



RNA viruses: The genomes are roughly 10kb, so

To get to 2x coverage, (some contigs): 200 reads at 100bp

to get 10x coverage (high confidence contigs): 1000 reads at 100bp.

100x coverage (major chunks of genome) 10000 reads at 100bp.

For Bacteria with 2Mb genome

To get 10x coverage at 100bp is 200K reads.

To get 100x coverage of course would take 2M reads,

Viruses:

Many of the HCV samples have a few K hits, so 10x suppression would allow you to sequence the genomes well (modulo scaffolding issues).

Many of them also had hepatitis B at a few tens to a few K reads.

HBV as a small genome of ~3kb so 10x suppression yields 30X + coverage.

Bacteria:

Samples with high counts: tens of K to low 100K.





10x suppression yields 20x coverage

Samples with Low counts, between 1500 and 10K counts.


10x suppression yields 1-5x coverage, yields some contigs, species info and possibly virulence islands.

100x suppression would be desirable, but you could do a lot with 20-50x suppression.

HAC normalization has been implemented in benchtop and automated systems



Suppression increases hits on pathogens in mock clinical samples



Multiple manifestations of the technique:


- Benchtop spin methods of normalization “HAC” and capture
- Manually operated capillary chromatography based methods
- Automated and multiplexed capillary based methods

Suppression enriches pathogen sequences and increases coverage




Suppression enriches for rare nucleic acids, including those derived from the pathogen


- Increases efficiency of SGS analysis
- Improves sequence coverage of pathogen




SGS “hits” on *F. tularensis* in human background (1:1000 spike-in)






Suppression increases miSeq hits on viral Pathogens in clinical samples: plasma



Increased coverage of viral genomes through suppression



Our bioinformatics pipeline enables efficient identification of pathogens in clinical samples



Nasopharyngeal swabs from patients
13 with respiratory complaints

SFAF June 7, 2012

The difference between short term & sustained productivity is the problem



**Short term areal production
of 30-50 gm/m²
–Commonly claimed**

**Long term areal production rates of
2 gm /m² (35 years of production)
–Ami Ben Amotz**




Sub optimal growth conditions: Irradiance, temp, salinity etc. and pond instabilities lead to lower long term production.

Closed PBRs can increase yield and postpone contamination, but have yet to be proved economically viable

Real time data on predator/pathogen load is enables proactive pond management

Presence of the biological agent can be necessary but not sufficient to crash




“Perhaps the most worrisome component of the large-scale algal cultivation enterprise is the fact that algal predators and pathogens are both pervasive and little understood.”

- DOE Draft Algal Biofuels Technology Roadmap (2009)

Parallels between this problem and that of the unknown bioterror were evident. We believed that elements of a RapTOR like approach would be effective.

Presumptive Identification of pond crash agents through sequencing



- ✓ **Probe independent: detection of unknowns**
- ✓ **Physical removal of host DNA prior to sequencing *up to 250 fold***
- ✓ **No dependence on microscopic analysis/expertise**
- ✓ **Effective for a wide variety of algal pathogens/predators**



| % hits | Best hit | type |
|--------|----------------------------------|-----------------------|
| 29 | <i>Cercomonas plasmodialis</i> ; | flagellate |
| 21 | <i>Aplanochytrium stocchinoi</i> | marine fungus |
| 16 | <i>Chaetonotus neptuni</i> | unsegmented worm-like |
| 4 | <i>Labyrinthuloides minuta</i> ; | plant pathogen |
| 4 | <i>Platyreta germanica</i> | parasite/predator |
| 4 | <i>Amphora cf. capitellata</i> | competitor |

SFAF June 7, 2012



Conclusions

We have developed molecular suppression techniques (HAC normalization, capture) that enable selective enrichment of pathogen NA in complex samples for targeted SGS analysis.

- **Suppression techniques implemented in benchtop (spin) and capillary chromatographic architecture**
- **in process of developing multiplexed automated capillary system**

We have demonstrated the effects of suppression on pathogen characterization in mock clinical and clinical samples

We have developed a metagenomics analysis pipeline that enables efficient identification of pathogen-derived NA sequences in complex clinical samples.



Acknowledgments

**Steve Branda
Ken Patel
Joseph Schoeniger**

**Michael Bartsch
Zachary Bent
Bryan Carson
Joshua Crawford
Deanna Curtis
Hanyoup Kim
Sadie LaBauvre
Pamela Lane
Stan Langevin
Chi-Chi May
Robert Meagher**

**Milind Misra
Nicholas Pattengale
Martha Perez-Arriaga
Amy Powell
Ron Renzi
Bryce Ricken
Anupama Sinha
Owen Solberg
Numrin Thaitrong
Mary Tran-Gyamfi
Victoria VanderNoot
Kelly Williams**