

# *Microbial forensics of select agents from trace samples; making the case for targeted sequencing*



5/29/2015 Tom Slezak, [slezak@llnl.gov](mailto:slezak@llnl.gov)

LLNL-PRES-672363

This work was performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under contract DE-AC52-07NA27344.

Lawrence Livermore National Security, LLC



# **Is there a role for targeted amplification NGS in Microbial Forensics?**

# **Is unbiased metagenomic sequencing the answer to every question?**

- **Clinical samples and environmental samples**
- **Detection vs. characterization vs. census vs. discovery vs. attribution...**

# **When can we draw a line between discovery of novel agents versus ID & characterization of known agents?**

- **Does every sequencing run have to address the potential of novel discovery?**
- **How does cost and time-to-answer factor in for the different use cases?**

**Targeted amplification gives you the precision of multiplexed PCR with the scalability of microarrays, plus the benefit of NGS as a readout providing ~100+bp of sequence/target depending on chemistry used.**

**Can target organism, gene, or SNP resolution, in any combination.**

**Up-front bioinformatics to target highly-informative regions at the desired resolution(s), instead of the random hits of unbiased metagenomic NGS.**

# What about host depletion and other clutter-mitigation techniques?

- These tend to be host-centric for clinical samples
- Not clear how this would work for environmental samples; it is all “clutter”!

***No magic bullet on the horizon for microbial forensics from environmental samples***

# Two basic types of targeted amplification:

- PCR-based
  - deeper amplification; hundreds to ~20K+ amplicons targeted
  - ~5hr prep
- Hybridization based
  - more highly scalable; more forgiving of sequence mis-matches; ~lower amplification
  - ~overnight hyb




# Many targeted amplification uses are possible:

- Symptomatic disease panels
- Biothreat agent panels
- AMR panel
- *Rapid forensic characterization panels*
- Etc.


# A pilot study showed us how well targeted amplification works

- ~100-1000x enrichment of the informative targeted regions
- Species-ID of a panel of Select Agents

Collaboration with Ken Frey of NMRC and Matt Dyer, Adam Allred, Brian Kelly of LifeTechnologies (now ThermoFisher)



# We observed 100-1000x enrichment of (informative) pathogen reads



**Can rapid initial forensic  
characterization be done from a  
trace sample using targeted  
amplification?**

# Consider a BioWatch agent panel:

- ~50 amplicons/oligos per agent to confirm genus and species
- ~50 amplicons/oligos per agent to query presence of key genes (virulence-related, gene-mediated resistance)
  - Additional amplicons/oligos to detect SNP-mediated resistance (e.g., Cipro, Doxy)
- ~500 amplicons/oligos per agent to establish SNP-based phylogenetic placement
  - 5x redundancy for ~100 SNPs to establish moderate-resolution phylogenetic placement.

*Potential to run directly from BioWatch filter extract...*

**Use of targeted amplification will permit sample multiplexing to greatly lower cost/sample of NGS characterization of known agents.**

***Obtain preliminary characterization and actionable information while you are doing your deep sequencing run or isolate culturing in parallel.***

***We need this in our biodefense toolkit!***