


Overview of the Requirements and Construction of Targets for Experiments on OMEGA and OMEGA EP









M. J. Bonino
University of Rochester
Laboratory for Laser Energetics

20th Target Fabrication Meeting Santa Fe, NM 20–24 May 2012

Summary

Complex targets have been developed for three experimental campaigns



- Cryogenic Compton radiography campaign
 - 60-beam OMEGA laser implodes a cryogenic capsule with 100- μ m-thick DT ice

Issue: tolerances on backlighter position relative to the capsule source

- Two-plamon-decay experiments
 - a 200-μm-diam molybdenum (Mo)-coated capsule co-centered with a GDP capsule

Issue: centering the two capsules to $\pm 50~\mu m$ (radius) and reducing the seam size


- Polar-drive shimmed (PDS) target
 - the PDS capsule has a constant inner radius and a varying wall thickness

Issue: wall-thickness variation and mounting strategy

Compton radiography cryogenic target assembly requires strict tolerances





 The capsule implosion is backlit with x rays (60 to 200 keV) to image the core Point-projection gold microwire positioning tolerances





Constructing the target assembly is a multi-step process













- 1. 250-μm-diam Be wire is formed at 400°C in a geometry avoiding beam interferences
- 2. $10-\mu$ m-diam gold wire array is assembled and characterized
- 3. Be wire, backlighter subassembly, and target holder are assembled on a fixture
- 14-μm-diam PBO fibers
 (Zylon) are constructed in a tripod arrangement

The capsule is positioned and joined to the target mount




• The completed assembly is measured using the Nikon software package





Cooling the Compton radiography target to 15 K induced dimensional changes (<100 μ m) that are manageable




This suggests a 200- μ m-diam core can be imaged.

Integrated two-plasmon-decay experiments on OMEGA will address electron coupling for direct-drive implosions




Symmetric or polar drive




- "Preheat" is from laser/ plasma interaction-generated electrons (50 to 150 keV) that can reduce compressibility
- The amount of preheat is measured by the K_{α} yield of a molybdenum capsule inside a 20- μ m-thick GDP capsule
- Requirements
 - ≥30 μm molybdenum on inner capsule
 - seam must be gap-free
 - concentricity to ±50-μm radius

Using a halfraum assembly approach, a molybdenum-coated capsule* is centered inside two hemispheres





1. Machine-matched sets of hemispheres (same wall thickness, measurement error $\pm 0.3~\mu m$; outer diameter, $\pm 5~\mu m$); to be completed at GA




2. Insert molybdenum-coated capsule using (a) formvar or (b) spider silk; add viscous glue to the seam to prevent laser light going into the center; to be completed at LLE

^{*}J. Jaquez, General Atomics.

A prototype target was made using formvar to center the inner capsule




Top view



218- μ m-diam glass capsule co-centered to 8 μ m of the GDP capsule outer diameter

Side view




Seam is in view, glass capsule is centered to 9 μ m of the outer diameter; seam width ~90 μ m

Mounting the inner sphere using spider silk is another viable approach to achieve ± 50 - μ m concentricity requirement




Top view



Glass capsule is within 20 μ m of center of the GDP capsule

Side view



Glass capsule is within 14 μ m of center of the GDP capsule


- Challenges: reduction of adhesive at seam
 - handling molybdenum coating on 140- μ m-diam glass




The fraction of hot electrons reaching the cold shell is measured using Mo-coated capsules



X-ray pinhole camera






These experiments currently provide an upper bound on the fraction of hot electrons that reach the cold shell (preheat).

Polar-drive shimmed (PDS) targets have been fielded on OMEGA






- 40-beam OMEGA shots
- Asymmetric intensity
- Displaced pointing


A shimmed capsule will implode uniformly, compensating for the asymmetric illumination.

A polar-drive shimmed target uses a capsule with variable wall thickness






Capsule thickness to compensate for 10% prolate intensity




Polar-drive shimmed targets are machining intensive






Transferring the capsule requires precision fixturing





Machined capsules are coated, characterized for gas retention, and fielded from a planar moving cryostat





Capsule sent from General Atomics, then Al sputter coated for gas retention



Advantages:

- Orientation
- Symmetry

Room-temperature moving cryostat with a spherical target (positioning accuracy <10 μ m, 1-h shot cycle)


Challenges: machining and accounting for orientation and gas retention

Results

Framed radiographs of polar-driven target implosions on OMEGA show improved implosion symmetry with target shimming



10 atm, D_2 -filled, 27- μ m-thick CH capsules, imploded with triple-picket pulses from 40 OMEGA beams, 14 kJ on target, backlit with Ti x-ray emission at 4.7 keV



The shimmed capsule implodes more uniformly with this beam pointing (120-, 140- μ m offsets for Rings 2 and 3).

Summary/Conclusions

Complex targets have been developed for three experimental campaigns



- Cryogenic Compton radiography campaign
 - 60-beam OMEGA laser implodes a cryogenic capsule with 100- μ m-thick DT ice

Issue: tolerances on backlighter position relative to the capsule source

- Two-plamon-decay experiments
 - a 200-μm-diam molybdenum (Mo)-coated capsule co-centered with a GDP capsule

Issue: centering the two capsules to $\pm 50~\mu m$ (radius) and reducing the seam size

- Polar-drive shimmed (PDS) target
 - the PDS capsule has a constant inner radius and a varying wall thickness

Issue: wall-thickness variation and mounting strategy