

Characterization of low-density HEDP laser target foam components using monochromatic soft x-ray sources

Martin Taccetti, Nick Lanier, Chris Hamilton

Los Alamos National Laboratory

Target Fabrication Meeting - Santa Fe, NM


May 22, 2012

We require accurate characterization of low-density foam components for successful radiation flow experiments on NIF and Omega

- Small deviations in foam density dramatically impact radiation propagation and complicate assessment of flow models
- Foam density is currently characterized in bulk and assumed to stay unchanged during machining
- Individual components are too light to weigh, thus normally do not make density measurements of single parts
- For this reason we constructed a monochromatic x-ray imager to measure the x-ray transmission of single component foams, used to derive a line-averaged density
- For best results, this measurement should be performed in conjunction with ICP-MS (inductively coupled plasma mass spectrometry) analysis, to constrain impurity levels
- Depending on the size and density of the part in question, the density measurement uncertainty can reach the limiting factor of the cold opacity uncertainty, $\sim 1\%$


We measure the transmission of soft x-rays through low-density foam components

- The x-ray imager consists of a soft x-ray generating device at either 2.3 or 5.4 keV, each with a dedicated sample chamber and CCD camera
- Recently re-configured to use contact radiography: X-rays are transmitted through the sample and imaged onto a CCD at $\sim 1\times$ magnification, recording line-averaged transmission as a function of position
- The soft x-rays are in the right regime to measure very small density variations in low-density C_8H_8 and SiO_2 foam targets
- We can measure line-averaged transmission to within $\sim 1\%$ relative uncertainty




Chromium $K\alpha$ (5.4 keV) setup

We illuminate a larger area of the detector by locating it outside the Rowland circle



The x-ray beam covers enough area on the detector to allow characterization of most HEDP target foams

2048 x 2048 pixel
CCD image



With knowledge of cold opacity, high resolution monochromatic transmission measurements constrain line-averaged density

Change in intensity as light passes through an absorptive medium

$$dI(\nu, z) = -\overset{\substack{\text{opacity} \\ \downarrow}}{k(\nu)} \overset{\substack{\text{intensity} \\ \downarrow}}{\rho(z)} I(\nu, z) \overset{\substack{\uparrow \\ \text{density}}}{dz} \overset{\substack{\uparrow \\ \text{distance}}}{dz}$$

Solving for intensity and integrating over path length

$$\begin{aligned} \int \frac{dI(\nu, z)}{I(\nu, z)} &= - \int k(\nu) \rho(z) dz \\ \ln I(\nu, z) - \ln I_0(\nu, z) &= - \int k(\nu) \rho(z) dz = -k(\nu) \bar{\rho} \int dz = -k(\nu) \bar{\rho} z \end{aligned}$$

line-averaged density
↓

thus $I(\nu, z) = I_0(\nu, z) e^{-k(\nu) \bar{\rho} z}$ (general formula)

We use a Dirac delta function for the monochromatic source and integrate over ν

$$\int I(\nu, z) \delta(\nu - \nu_0) d\nu = \int I_0(\nu, z) e^{-k(\nu) \bar{\rho} z} \delta(\nu - \nu_0) d\nu$$

or

$$I(\nu_0, z) = I_0(\nu_0, z) e^{-k(\nu_0) \bar{\rho} z}$$

mono-energetic
cold opacity

And then solve for the line-averaged density

$$-\ln \frac{I(\nu_0, z)}{I_0(\nu_0, z)} = -\ln T(\nu_0, z) = k(\nu_0) \bar{\rho} z \quad \therefore \quad \bar{\rho} = -\frac{\ln T(\nu_0, z)}{k(\nu_0) z}$$

Taking partial derivatives obtain error terms for T, k, z

Line-averaged density $\bar{\rho} = -\frac{\ln T(\nu_0, z)}{k(\nu_0)z}$

Take partial derivative to obtain each error term

$$\begin{aligned}\partial \bar{\rho} &= \left| \frac{\partial \bar{\rho}}{\partial T} \right| \partial T + \left| \frac{\partial \bar{\rho}}{\partial k} \right| \partial k + \left| \frac{\partial \bar{\rho}}{\partial z} \right| \partial z \\ &= \left| \frac{1}{T(\nu_0, z) k(\nu_0) z} \right| \partial T + \left| \frac{\ln T(\nu_0, z)}{k(\nu_0)^2 z} \right| \partial k + \left| \frac{\ln T(\nu_0, z)}{k(\nu_0) z^2} \right| \partial z\end{aligned}$$

Transmission
error term

Opacity
error term

Axial Distance
error term

For opacity and axial distance we can assume the following uncertainties:


Cold opacity $\rightarrow \frac{\partial k}{k(\nu_0)} = 1\%$ and Distance $\rightarrow \partial z = 1 \mu\text{m}$

The experimental error in the x-ray transmission measurement is dominated by photon statistics

Sample
Specifications

SiO_2
 $\rho \sim 125 \text{ mg/cc}$
 $h = 200 \mu\text{m}$
 $d = 2000 \mu\text{m}$


$m \sim 78 \mu\text{g}$



All subsequent results are extracted from
250 accumulated images

The pixel-to-pixel uncertainty is $\sim 1\%$, greater reductions are achieved by averaging over larger regions of interest

Pixel-to-pixel variation for 250 frame ensemble



A $75\text{-}\mu\text{m}$ by $75\text{-}\mu\text{m}$ region of interest yields a $\sigma_T \sim 0.2\%$


Densities inferred with our x-ray assembly are statistically equivalent to measurements made at Brookhaven's NSLS

	X-ray Assembly	NSLS	~ 15 μm Silicon Standard
Transmission @ 5.415 keV	$59.4 \pm 0.8 \%$	$59.0 \pm 0.5 \%$	

Sample Specifications

SiO_2
 $\rho \sim 70 \text{ mg/cc}$
 $h = 200 \mu\text{m}$
 $d = 400 \mu\text{m}$


$m \sim 2 \mu\text{g}$



A blind study provided the best option to challenge our foam manufacture and characterization techniques

- Multiple SiO₂ samples at three desired densities were manufactured
- Each sample was 200 microns thick with a 400 micron diameter
- Densities would range between 10-160 mg/cc and were cast in place
- Since overall part mass was too small for gravimetric measurements, bulk witness samples were manufactured and characterized
- Samples were also characterized on the x-ray DCS
- A subset of samples would also be characterized at Brookhaven National Laboratory's NSLS


Blind study results show the bulk density often differs significantly from individual samples



- Bulk density measurements for witness samples disagree with that derived from x-ray transmission for small cast in place SiO_2 foams
- Significant variability existed between samples of supposedly identical density
- The x-ray DCS clearly resolves this variability

Increased water content can lead to higher inferred densities if not properly accounted for


X-ray assembly data showing the hygroscopic nature of SiO_2 foams





After baking, the densities approach their original values.

Foam density characterization station is already impacting selection of NIF HEDP targets

- Our x-ray assembly is impacting the first set of non-ignition related experiments being performed at NIF in support of development and validation of models used in our codes
- HEDP targets include 125 mg/cc silica aerogel or chlorinated plastic foams that require detailed characterization



Some targets, such as this HiPE ($\text{C}_8\text{H}_7\text{Cl}$) one - were eliminated from our shot sequence after characterization



Averaging over ROI:

$$\rho_{\text{avg}} = 124.0 \pm 1.3 \text{ mg/cm}^3$$

$$(\rho_{\text{witness}} = 126 \text{ mg/cm}^3)$$




Assumptions:

HiPE cold opacity @5415 eV
 $= 90.3 \text{ cm}^2/\text{g} \pm 1\%$

thickness error is $\pm 1 \mu\text{m}$

The x-ray assembly also provides a good test of density uniformity

$$\rho_{\text{avg}} = 124.0 \pm 1.3 \text{ mg/cm}^3$$



Summary

- We have constructed a monochromatic x-ray imager to quantify the density of low-mass single component foams
- Density inferences from the x-ray assembly are identical to those obtained at NSLS
- Small cast-in-place SiO_2 foams often exhibit high variability in their density, and are poorly represented by bulk measurements on witness samples
- For the masses and sizes tested, the x-ray assembly resolved line-average transmission down to $\sim 1\%$
- When left exposed to air, the SiO_2 samples absorbed water, leading to an increase in inferred densities
 - easily resolved by the x-ray DCS
 - more characterization of this effect is essential to assure foam components at the time of experiments are understood
- Current efforts implementing 2.3 and 8.0 keV sources are underway