

Ultrasensitive Detection Using an Optical Trap

D. J. Vieira, M.D. Di Rosa, X. Zhao
S.G. Crane, M. Hausmann, and J. Wu

- Near resonant absorption (high selectivity)
- 10^7 “scattered” photons / sec / atom (high sensitivity)
 - single atoms have been detected in a MOT

Application: trace isotope determination is important for environmental and nonproliferation / treaty verification
- consider Cs isotopes:

Cs133 100%	Cs134 2.1 a □.605, ...	Cs135 2.3x10 ⁶ a no □	Cs136 13 d □.818, ...	Cs137 30 a □.117, ...	
Xe132 27%	Xe133 5.2 d □.081, ...	Xe134 10%	Xe135 9.1 h □.250, ... □.2.6 x 10 ⁻⁶	Xe136 8.9%	Xe137 3.8 m □.455, ...
		6.7%	6.5%	6.2%	

Fission Yield


Large $^{135}\text{Xe}(n, \gamma)$ cross section diverts yield to A=136 mass chain

- $^{135}\text{Cs}/^{137}\text{Cs}$ isotopic ratio in reactor fuel is sensitive to neutron fluence history

Optical traps have the potential of achieving much greater sensitivities

- 100x better than existing techniques (TIMS, AMS, RIMS) depending on ΔZ

Ultra-Sensitive Detection Using an Optical Trap Coupled to a Mass Separator






Figure 9

Di Rosa *et al.* Appl. Phys. B (in press)



Summary & Outlook

- We have demonstrated the measurement of $^{135}\text{Cs}/^{137}\text{Cs}$ at ultrasensitive detection levels.
 - 4,000 atom detected / overall efficiency = 0.5%
 - 1 million atoms per sample
 - linear over 3 orders of magnitude with an isotope ratio accuracy of better than 10%
 - high selectivity achieved by combining mass separator & optical trap ($>10^{12}$ suppress of ^{133}Cs relative to $^{135,137}\text{Cs}$)
- This method can be future improved & extended to other species
 - reduce implantation foil “memory”
 - increase detection sensitivity to the few atom level
 - Sr and Kr radioactive isotopes
- Real world applications
 - nonproliferation, threat reduction and the environment

