





# The split-GFP system

# Create a new detection system based on protein complementation




# Finding complementing pairs of GFP



Only (GFP 10-11 + GFP 1-9) and (GFP 11 + GFP 1-10) pairs from Superfolder GFP showed complementation and fluorescence.


# The Split-GFP tagging system



*Cabantous S. et al. (2005). "Protein tagging and detection using engineered self-assembling fragments of green fluorescent protein". Nat. Biotechnology 23(1), 102-7.*

# Improving GFP 1-10 solubility

1



2

GFP 1-10 LIB  
(pET T7)

Expression strain  
SR-GFP 11 (*tet*)

3


Induction  
IPTG,  
Tetracycline

(1) Mutate Gene

Recombine  
Optima

(2) Clone

(3) Select  
Soluble clones




***In vivo* complementation  
of the two GFP fragments**





# The large fragment 1-10 was engineered to improve performance and solubility

• *In vitro*



• *In vivo*



Small GFP 11 tag was engineered to minimize its effect on passenger folding and solubility



# *In vitro* complementation

Express X-S11

Fractionate  
& Centrifuge


+ 9M Urea

soluble

insoluble ©

+ S1-10

15  
min





**A calibration © with a standard protein enables the quantification of soluble protein from assay fluorescence**



# *In vivo complementation*

2 independently-inducible plasmids in *E. coli*



*Co-induction of tagged protein & detection fragment*  
*Sequential induction of tagged protein then detection fragment*

*total expression* →  
*soluble* →

# Applications of Split-GFP

- High throughput Expression Screening: soluble, insoluble, total expression
- Assay of proteins during workup & purification
- Refolding assays: find optimal refolding conditions
- Find soluble mutants of a protein (directed evolution)
- Find soluble domains of a protein