

Raman spectroscopy calibration/analysis/errors – large scale containers

John M. Berg
LANL, NMT-11
October 22, 2002

Topics


- **Brief overview of the measurement method.**
- **Measurement results to date on large cans.**
- **Methods used to convert raw data to partial pressures.**
- **Detection thresholds and how they are defined and determined.**

Overview of Raman application

- Raman spectroscopy scatters light off of a sample and measures the scattered intensity as a function of wavelength.
- Each molecule scatters at characteristic wavelengths which are well known, allowing straightforward identification.
- Intensity is proportional to concentration.
- Fiber optics deliver and collect light.

Calibration gas mixture in Raman sample chamber

10 min. at 0.5 watts




Composition	Detection Threshold (torr)
5.02 torr CO_2	0.5
5.00 torr CO	0.9
7.02 torr H_2	0.5
4.99 torr CH_4	0.15
3.49 torr O_2	0.5
64.5 torr N_2	1.2
4.98 torr Ar	N/A
4.99 torr He	N/A


Methods of data treatment

- Raw data contain artifacts due to ionizing radiation hitting the detector (cosmic rays and local background). These are rare events that are removed by threshold filtering.
- Noise on inevitable background signal is reduced by signal averaging. It is the limiting factor in determining sensitivity.
- Detection thresholds are defined as the gas pressure that would produce a peak 5 times the amplitude of the background signal standard deviation.

Raman data treatment




Raw spectrum, 1 minute



Average of 10, outliers discarded


Detector background subtracted


Typical Raman data



- Calibration gas included H_2 , O_2 , N_2 , CO and CO_2 .
- Gas signals do not vary appreciably between probes.

Container 1 shows no Raman-detectable gases after > 9 months.






Raman Detection Thresholds for Container 1

Gas pressure that would produce a peak 5 times the amplitude of the background signal standard deviation in 10 minutes.

	Torr	psi	std. atm.
H₂	3	0.06	0.4 %
O₂	4	0.08	0.5 %
CO₂	4	0.08	0.5 %
N₂O	4	0.08	0.5 %
N₂	10	0.08	1.3 %

Calibration Gas Mixture in Empty Can

- Normalized to N_2 known value of 4.03%.
- Standard deviations of 5 to 8 %.



Application to Containers

- Raman is in-situ measurement, consuming no gas.
- Measurement frequency not limited by need to conserve sample.
- Initial measurements taken several per day.
- Measurement rate tapering to once per month after 3 months.

Summary

- Pure oxide in Container 1: no Raman signals from detectable gases (essentially everything but noble gases).
- Detection thresholds for 10 minute acquisitions are conservatively estimated at 3 torr for H_2 , 4 torr for O_2 , CO_2 and N_2O , and 10 torr for N_2 .