SFE efforts in FY02 and status of technical basis

W. Kirk Hollis and Kelly Hakonson







Status

- Re-run of selective MIS samples Complete
 - As Received
 - 24 hour reabsorption
 - 7 day reabsorption
- Report in final stages
 - SFE, IGA and selective TGA results
 - Sample protocols used
 - Need to compare SFE with TGA for 6 samples







Sample Categories

- Same as Gillispie report
 - Pure PuO2
 - Pu/U mixed
 - Mg and Ca
 - Na and K
 - Others Impurities
- 60+ samples run







Definitions

- False Positive: IGA > SFE and exceeded 0.5%
- TGA results are reported as LOI up to 500°C and small sample size (mgs)
- 25 Relative Percent Difference (RPD) between SFE and IGA is consider acceptable







- Pure PuO₂
 - No False Positives
 - 11 samples run
 - 8 SFE/IGA <25 RPD
 - 2 IGA > SFE {MISNE2 (0.08, 0.03 @62.5%) and B1579 (0.17, 0.09@47.1%) }
 - BLO39 SFE <0.1% , IGA 0.41-0.44, TGA 0.51-0.57
 - Review of TGA results needed for final determination







- Pu/U mixtures
 - No False Positives
 - 10 samples have been run
 - 6 SFE/IGA <25 RPD
 - 2 IGA > SFE at 24 hours, same samples are part of 6 above at 7 days
 - SPC-711 (non calcined) IGA X2 SFE (0.4 / 0.2)
 - SFE and IGA agreed following calcination (<0.03%)
 - TGA results under review







- Mg and Ca
 - 1 false positive (7 day, J2201, 0.35% and 0.50%)
 - Suggestion in lower moisture control value to 0.25% to address this concern
 - 26 samples run
 - 14 SFE/IGA <25RPD
 - 13 SFE/IGA >25RPD
 - 5 moisture was >1% and observation reported heterogeneity
 - 2 were part of multiple analysis and average of the 3 run meet criteria
 - 4 samples form Hanford (R437 and 440) initial results showed agreement with SFE/IGA and TGA. Later results > criteria.
 Reviewing TGA results for better understanding.







- Na and Ca
 - No False Positives
 - 7 samples run
 - 4 SFE/IGA <25RPD
 - 3 SFE/IGA >25RPD (all SFE higher)
 - One would be considered False Positive for IGA (24 hour)
 - Run log reported plugging of restrictor
 - Restrictor replaced (7 day run of same sample passed)







- Other Impurities
 - No False Positives
 - 9 samples run
 - 5 SFE/IGA <25RPD
 - 4 SFE/IGA >25RPD
 - 2 results of both <0.1%</p>
 - 2 SFE > IGA
 - » 1 IGA false positive (0.65% and 0.44%)
 - » Both samples run log (inhomogeneous)







Current Status

- TGA input to finalize report for review process
- Instrument has not been used for 6-9 months and downtime has cause operational issues, currently undergoing repair.
- New SFE instrument should be installed in PF-4 sometime in this FY





