

Materials Identification and Surveillance – Gas Chromatography Analysis of 3013 Storage Cans 1 and 10

David Harradine, Rhonda McInroy (C-PCS)

Kirk Hollis (C-ACT)


John Berg, Max Martinez, Jim McFarlan, Kirk Veirs,
Laura Worl (NMT-11)

Dallas Hill, Jim Stewart (ESA-EPA)


94-1 R&D End of Year Technical Review
October 21 & 22, 2002

GC Sampling Procedures & Sources of Error

- GC sampling completely automated using LabView software.
- Expansion factor from sampling volume to GC ~ 40 . A limitation to sensitivity.
- Absolute GC measurement turned into a relative percentage and multiplied by measured can pressure for partial pressure of various gases.




Agilent 6890 Configuration & Sources of Error cont.




- Thermoconductivity (TCD) detector.
- Molsieve & Plot-Q columns in series.
- Argon used as carrier gas for H_2 sensitivity.
- Sensitivity to other gases decreases.
- Peak overlap between He and H_2 .

Typical Gas Chromatograph, Agilent 6890



Error in Measurement and Detection Limits




Analysis error arises from:

- Measurement error in pressure
- Run to run variations in GC
- Integration errors
- Calibration errors

QA controls:


- 6-12 month multipoint calibration of GC
- single point calibration for each day of use
- Use of “PeakFit for data analysis
- Test of sampling procedure using Can 10

Species	Error (\pm Torr)	Detection Limit (Torr)
CO ₂	3.8	2.0
H ₂	1.3	0.5
O ₂	2.3	2.0
N ₂	2.8	2.0
CH ₄	1.2	0.7
CO	4.7	2.0



PeakFit Integration

Data analysis involves a multiparameter fit to the baseline followed by curve fitting using either a Haarhoff-Van der Linde or an Exponentially Modified Gaussian function.



GC Monitoring 3013 Storage Can 1


Can 1, sintered at 950°C and back filled with Helium.

Pressure lowered at 50 days.

Data are corrected for



- temperature (average of can thermocouples).
- zero point drift of pressure gauge.

GC analysis is relative and multiplied by measured can pressure.




GC Monitoring 3013 Storage Can 1

Days	CO ₂	N ₂ O	He	H ₂	O ₂	N ₂	CH ₄	CO
0.0	0	0	19.9	0	0	0	0	0
0.1	0	0	20.1	0	0	0	0	0
3.9	0	0	20.1	0	0	0	0	0
6.9	0	0	20.1	0	0	0	0	0
21.0	0	0	19.9	0	0	0	0	0
34.8	0	0	19.3	0	0	0	0	0
48.9	0	0	12.9	0	0	0	0	0
62.9	0	0	12.5	0	0	0	0	0
77.0	0	0	12.7	0	0	0	0	0
250.9	0	0	12.5	0	0	0	0	0



Conclusion:
Nothing of consequence
in happening in Can 1

GC Monitoring 3013 Storage Can 10



Purpose of Can 10:

- Reactivity of gases with walls.
- Sample procedure check.
- Pressure rise between 100 and 350 days may be real.
- Change in composition is close to our limits of detection.

Differences in % Composition from Expected

Species	Volume %
CO ₂	3.98%
H ₂	3.95%
O ₂	4.01%
N ₂	4.03%
CH ₄	4.08%
CO	4.07%
Ar	75.88%

Argon is used as the carrier gas and can't be measured, assumed constant.

