Materials Identification and Surveillance Small Scale Reactors

David Harradine, Rhonda McInroy (C-PCS) Kirk Hollis (C-ACT)

John Berg, Max Martinez, John Morris, Kirk Veirs, Laura Worl (NMT-11) Dallas Hill, Jim Stewart (ESA-EPA)

3013 Monitoring and Surveillance Quarterly Meeting February 25-27, 2003\
Savannah River Site





Small Scale Reactors Outline


- •Installation in glove box and Work Authorization
- •The PF-4 Noise Problem
- •Initial sample selection
- •Results
- •Why results are different than Duffey and Livingston
- •Future work






Small scale container surveillance samples

- Containers are in heated blocks (45 locations) with heated gas manifold
- ♦ Matrix assembly installed in plutonium glove box, Sept. 2002
- → Plumbing complete, facility leak-tested, Dec. 2002
- First sample introduced and equipment checkout commenced, Dec. 2002
- → Grandé Noise discovered, modifications tested, Jan. 2003
- ❖ Process Hazard Analysis modified, Work Instruction completed Feb. 2003
- ♦ Work and Worker Authorization and initial samples introduced Feb. 13, 2003


















Grandé Noise in GC spectrum Problem diagnosed in a few days, two weeks to implement fix







Initial Sample Selection Goals

- ♦ Compare results from first Large Container material to small container configuration. Use same oxide material, keep dry, and heat to 55°C. (But did not calcine material again-material spent 12 months exposed to glove box atmosphere)
- ♦ Measure reproducibility of entire approach. Have five samples exactly the same. Chose to to have these samples at 0.5 wt% water.
- Compare to Duffey and Livingston's work. One sample at 2 wt% water.





Material Characteristics

PEOF1

87.81% Pu

0.09% LOI (Feb 7, 2002)

 $1.08 \text{ m}^2/\text{g SSA}$


11.5 g cm⁻³ pycnometer density

location	mass	water added	fill gas
C3	9.7 g	194 μΙ Η ₂ Ο	He
C4	10.2 g	52 μl H ₂ O	He
C5	9.4 g	47 μl H ₂ O	He
C6	9.1 g	45.5 μl H ₂ O	He
C7	10.2 g	52 μl H ₂ O	He
C8	9.3 g	47 μl H ₂ O	He
C9	9.6 g	0 μl H ₂ O	He
		-	







Calibration of gases with high GC sensitivity




Calibration of gases with low GC sensitivity




GC results show significant changes in gas composition



Temperature in Small Scale Containers is Stable



Pressure in Small Scale Containers



Calculation Methodology Using Container C3 as an Example

 CO_2 N_2O He H_2 O_2 N_2 CO

Integrated Area				
Feb. 13 Feb. 18				
0	1769			
0	1486			
259705	325074			
0	3467			
25	25			
117	809			
0	200			

1					
F	Р				
(To	orr)				
Feb. 13	Feb. 18				
0.00	1.06				
0.00	0.89				
19.75	24.72				
0.00	0.15				
0.01	0.01				
0.05	0.38				
0.00	0.10				

Moles in container				
Feb. 13	Feb. 18			
0	5.53E-06			
0	4.65E-06			
0.000141	0.00013			
0	7.85E-07			
7.27E-08	5.35E-08			
3.88E-07	1.97E-06			
0	4.99E-07			

Change		
in moles		
10 ⁻⁶ moles		
5.53		
4.65		
0.79		
-0.02		
1.58		
0.50		

Rate of
Change
nmol s ⁻¹ W ⁻¹
0.622
0.522
0.088
-0.002
0.178
0.056

Sum 19.81 27.30 Gauge 19.18 27.57 Sum of O 16.2 Sum of N 12.5 Sum of C 6.0

Moles of gas in one monolayer 132





Results of Five Days Surveillance

rate of change in nanomoles s ⁻¹ W ⁻¹						
Location	CO ₂	N_2O	H ₂	O_2	N_2	CO
C3 C4 C5 C6	0.622	0.522	0.088	-0.002	0.178	0.056
C4	0.619	0.647	0.098	-0.007	-0.076	0.074
C5	0.702	0.619	0.096	-0.004	0.174	0.073
C6	0.591	0.645	0.094	-0.007	0.162	0.031
C7	0.802	0.646	0.159	-0.001	0.189	0.023
C8	0.557	0.593	0.078	-0.185	-0.348	0.042
C9	0.195	0.080	0.012	0.000	0.196	0.050

Average	0.678	0.639	0.112	-0.005	0.112	0.050
Sigma	0.095	0.013	0.032	0.002	0.126	0.027





Comparison of Initial Results to Duffey and Livingston

H ₂ generation rate			
	Duffey and Livingston	This Work	
	nmol W ⁻¹ s ⁻¹		
2 wt%	2.1 ^a	0.088 ^b	
1.0 wt%	0.056 ^a		
0.5 wt%	0.039 ^b	0.112 ^b	
0.0 wt%		0.012 ^b	

^a calcined at 700C

Duffey and Livingston experiments at room temperature.

This work at 55C with significant temperature gradients. Liquid water likely formed on the lid, which was cooler than the material. Therefore, amount of water adsorbed onto the material is governed by the vapor pressure of water at the coolest spot, which is likely the same for 0.5wt% and 2.0wt% experiments. Thermal modeling of 3013 containers by Hensel and by Knight indicate that large temperature gradients will exist in real containers and the container walls will be the coolest areas.





^b calcined at 950C

Future Work

- •Plutonium oxides of different surface areas with varying amounts of water, additional duplicates.
- •Salts with varying amounts of water, duplicates.
- •Laura's talk and discussion.





Small Scale Surveillance Experimental Design

Container Features

- 5 ml internal volume
- Scaled to 1/500th of 3013 Storage Can
- Heated to appropriate T seen in full scale containers.

Material

- Locations contain MIS materials that meet the 3013 packaging criteria and some materials that do not meet the 3013 criteria (>0.5 wt.% moisture, calcination T)
- Initial 6 months include 5 sample blanks and 8 sample duplications
- Variables include moisture content and form of moisture (hydrated salts, hydroxides, physisorbed water), fill gas, temperature



FY2002 Modifications

- ⇒ Space and safety requirements led to container lid modifications
 - More robust design for operation inside glovebox (vs. open front hood) that include new valves, pressure transducer
 - Contamination release minimized
- ⇒ GC Modifications provided increase in sample sensitivity by factor of 100.





Small Scale P&ID Diagram GCV-1 45 μl sampling volume Gas Chromatograph HV-A1a SHV-A2a SHV-A3a SHV-A4a SHV-A5a SHV-A6a SHV-A7a SHV-A8a SHV-A9a Ashcroft 100 pressure Optional external GC psia pressure gauge sample introduction transducer Manifold Bo X Four port flow selector Position 2 ⊗HV-B2b ⊗HV-B3b ⊗HV-B4b ⊗HV-B5b ⊗HV-B6b ⊗HV-B7b ⊗HV-B8b ⊗HV-B9b B1 Vacuum Manifold Co (X) HV-C1a ⊗HV-C2a ⊗HV-C3a ⊗HV-C4a ⊗HV-C5a ⊗HV-C6a ⊗HV-C7a ⊗HV-C8a ⊗HV-C9a 10 Multiposition Dead-end flowpath SD configuration Valco Valve HV-C1b ⊗HV-C2b ⊗HV-C3b ⊗HV-C4b ⊗HV-C5b ⊗HV-C6b ⊗HV-C7b ⊗HV-C8b ⊗HV-C9b Optional internal GC =Manifold D: X Manifold E-0 X 00 HV-E1a ⊗HV-E2a ⊗HV-E3a ⊗HV-E4a ⊗HV-E5a ⊗HV-E6a ⊗HV-E7a ⊗HV-E8a ⊗HV-E9a ⊗HV-E4b ⊗HV-E5b ⊗HV-E6b ⊗HV-E7b ⊗HV-E8b ⊗HV-E9b Calibration maintained to assure calibration on the 45 reactor Low dead space manual valves gauges.

Los Alamos